

c-Kit Rabbit pAb


ABZ1695

Product Overview


Name	c-Kit Rabbit pAb
Accession(Primary)	P10721
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 110 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (P10721)


Gene Name	KIT
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded Human Mammary cancer using c-Kit antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using c-Kit antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human Mammary cancer using c-Kit antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.

