

14-3-3 alpha/beta Rabbit mAb

ABZ2569

Product Overview

| | |
|-----------------------------------|-------------------------------|
| Name | 14-3-3 alpha/beta Rabbit mAb |
| Catalog # | ABZ2569 |
| Clonality | Monoclonal |
| Accession(Primary) | P31946 |
| Application Note (Approx.) | WB1:500-1:1000 IHC0.736111111 |
| Precautions | |


Target information(P31946)

Synonyms**Gene ID** 7529**Other Names****Function****Cellular location****Note**



Western blot analysis of 14-3-3 alpha/beta in rat brain, rat lung lysates using 14-3-3 alpha/beta antibody.

Western blot analysis of 14-3-3 alpha/beta in HeLa, A549, MCF-7, U251, C2C12 lysates using 14-3-3 alpha/beta antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using 14-3-3 alpha/beta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

