

# A1BG Rabbit pAb

ABZ1668

## Product Overview

<b>Name</b>	A1BG Rabbit pAb
<b>Catalog #</b>	ABZ1668
<b>Clonality</b>	Polyclonal
<b>Accession(Primary)</b>	P04217
<b>Application Note (Approx.)</b>	WB1:500-1:2000 IHC1:100-1:300
<b>Precautions</b>	

## Target information(P04217)

**Synonyms****Gene ID** 1**Other Names****Function****Cellular location****Note**



Western blot analysis of A1BG in K562 lysates using  
A1BG antibody.



Immunohistochemistry analysis of paraffin-embedded Human ovary using A1BG antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human ovary using A1BG antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human uterus using A1BG antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human uterus using A1BG antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

---