

# ABCG2 Rabbit pAb


ABZ1281

## Product Overview

<b>Name</b>	ABCG2 Rabbit pAb
<b>Catalog #</b>	ABZ1281
<b>Clonality</b>	Polyclonal
<b>Accession(Primary)</b>	Q9UNQ0
<b>Application Note (Approx.)</b>	WB1:500-1:2000 IHC1:100-1:300
<b>Precautions</b>	

## Target information(Q9UNQ0)


**Synonyms****Gene ID** 9429**Other Names****Function****Cellular location****Note**



Western blot analysis of ABCG2 in various lysates using ABCG2 antibody.



Western blot analysis of ABCG2 in various lysates using ABCG2 antibody.



Western blot analysis of ABCG2 in HT-29 lysates using ABCG2 antibody



Immunohistochemistry analysis of paraffin-embedded Human brain using ABCG2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of ABCG2 in HT-29 lysates using ABCG2 antibody.