

KD-Validated Anti-ABCB1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1002

Product Overview

Name	KD-Validated Anti-ABCB1 Rabbit Monoclonal Antibody
Accession(Primary)	P08183
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human P Glycoprotein 1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 141 kDa , observed, 130-180 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:00-1:1,000

Target information (P08183)

Gene Name	ABCB1
About Gene Symbol	ABCB1 (ATP Binding Cassette Subfamily B Member 1) is a cell surface marker involved in adaptive immunity. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

