

KD-Validated Anti-TMEM173 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1003

Product Overview


Name	KD-Validated Anti-TMEM173 Rabbit Monoclonal Antibody
Accession(Primary)	Q86WV6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TMEM173
Species Reactivity	Human
Calculated MW	Predicted, 42 kDa, observed, 35 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (Q86WV6)


Gene Name	STING1
About Gene Symbol	STING1 (Stimulator Of Interferon Response CGAMP Interactor 1) is an innate immune sensor involved in JAK-STAT/cytokine signaling and innate immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely




Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HT-1080 cells with TMEM173 antibody (Cat#ABG1003, 1:1,000). Nuclei were stained blue with DAPI; TMEM173 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.