

KD-Validated Anti-Peptidylprolyl isomerase B Rabbit Monoclonal Antibody


ABG1006

Product Overview


Name	KD-Validated Anti-Peptidylprolyl isomerase B Rabbit Monoclonal Antibody
Accession(Primary)	P23284
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Cyclophilin B
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 24 kDa , observed, 18 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P23284)


Gene Name	PPIB
About Gene Symbol	PPIB (Peptidylprolyl Isomerase B) is a peptidylprolyl isomerase (cyclophilin) involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC