

KD-Validated Anti-DSTN Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1008


Product Overview

Name	KD-Validated Anti-DSTN Rabbit Monoclonal Antibody
Accession(Primary)	P60981
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Destrin
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 19 kDa, observed, 17 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P60981)

Gene Name	DSTN
About Gene Symbol	DSTN (Destrin, Actin Depolymerizing Factor) is an actin regulatory protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of H9C2 cells with anti-DSTN antibody (Cat#ABG1008, 1:1,000).

Nuclei were stained blue with DAPI; DSTN was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.