

KD-Validated Anti-Laminin gamma 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1009

Product Overview


Name	KD-Validated Anti-Laminin gamma 1 Rabbit Monoclonal Antibody
Accession(Primary)	P11047
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Laminin gamma 1
Species Reactivity	Human
Calculated MW	Predicted, 178 kDa, observed, 220-250 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P11047)


Gene Name	LAMC1
About Gene Symbol	LAMC1 (Laminin Subunit Gamma 1) is an extracellular matrix laminin protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC