

KD-Validated Anti-Histone H1.0 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1012

Product Overview


Name	KD-Validated Anti-Histone H1.0 Rabbit Monoclonal Antibody
Accession(Primary)	P07305
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Histone H1.0
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 21 kDa , observed, 28 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P07305)


Gene Name	H1-0
About Gene Symbol	H1-0 (H1.0 Linker Histone) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genulin Biotechnologies LLC



Copyright ©2024 Genulin Biotechnologies LLC



Copyright ©2024 Genulin Biotechnologies LLC