

KD-Validated Anti-IQGAP2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1021

Product Overview

Name	KD-Validated Anti-IQGAP2 Rabbit Monoclonal Antibody
Accession(Primary)	Q13576
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human IQGAP2
Species Reactivity	Human
Calculated MW	Predicted, 181 kDa , observed, 181 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (Q13576)

Gene Name	IQGAP2
About Gene Symbol	IQGAP2 (IQ Motif Containing GTPase Activating Protein 2) is a GTPase-activating protein involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

