

KD-Validated Anti-Dynamin 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1032

Product Overview


Name	KD-Validated Anti-Dynamin 1 Rabbit Monoclonal Antibody
Accession(Primary)	Q05193
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Dynamin 1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 97 kDa, observed, 97 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (Q05193)


Gene Name	DNM1
About Gene Symbol	DNM1 (Dynamin 1) is a dynamin GTPase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blotting analysis using anti-dynamin 1 antibody (Cat#ABG1032). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-dynamin 1 antibody (Cat#ABG1032, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-dynamin 1 antibody (Cat#ABG1032). Dynamin 1 expression in wild type (WT) and dynamin 1 (DNM1) shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-dynamin 1 antibody (Cat#ABG1032, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Immunocytochemical staining of C2C12 cells with anti-dynamin 1 antibody (Cat#ABG1032, 1:1000). Nuclei were stained blue with DAPI; Dynamin 1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 µm.