

KD-Validated Anti-ITPK1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1042

Product Overview

Name	KD-Validated Anti-ITPK1 Rabbit Monoclonal Antibody
Accession(Primary)	Q13572
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human ITPK1
Species Reactivity	Human
Calculated MW	Predicted, 46 kDa, observed, 42 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000


Target information (Q13572)

Gene Name	ITPK1
About Gene Symbol	ITPK1 (Inositol-Tetrakisphosphate 1-Kinase) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-ITPK1 antibody (Cat#ABG1042). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-ITPK1 antibody (Cat#ABG1042, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-ITPK1 antibody (Cat#ABG1042). ITPK1 expression in wild type (WT) and ITPK1 shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-ITPK1 antibody (Cat#ABG1042, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of ITPK1 expression in HepG2 cells using ITPK1 antibody (Cat#ABG1042, 1:2,000). Green, isotype control; red, ITPK1.



Immunocytochemical staining of HepG2 cells with anti-ITPK1 antibody (Cat#ABG1042, 1:1,000).

Nuclei were stained blue with DAPI; ITPK1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.