

KD-Validated Anti-MERTK Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1046


Product Overview

Name	KD-Validated Anti-MERTK Rabbit Monoclonal Antibody
Accession(Primary)	Q12866
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human MERTK
Species Reactivity	Human
Calculated MW	Predicted, 110 kDa , observed, 188 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (Q12866)

Gene Name	MERTK
About Gene Symbol	MERTK (MER Proto-Oncogene, Tyrosine Kinase) is a receptor tyrosine kinase involved in receptor-mediated signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HAP-1 cells with MERTK antibody (Cat#ABG1046, 1:1,000). Nuclei were stained blue with DAPI; MERTK was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.