

KD-Validated Anti-Inhibin subunit beta B Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1062

Product Overview


Name	KD-Validated Anti-Inhibin subunit beta B Rabbit Monoclonal Antibody
Accession(Primary)	P09529
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Inhibin beta B
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 45 kDa, observed, 52 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P09529)


Gene Name	INHBB
About Gene Symbol	INHBB (Inhibin Subunit Beta B) is a cytokine and immune signaling protein involved in TGF- β signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HT-1080 cells with anti-inhibin subunit beta B antibody (Cat#ABG1062, 1:1,000). Nuclei were stained blue with DAPI; Inhibin subunit beta B was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 μ m.