

KD-Validated Anti-CRABP1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1073

Product Overview

Name	KD-Validated Anti-CRABP1 Rabbit Monoclonal Antibody
Accession(Primary)	P29762
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CRABP1
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 16 kDa, observed, 14 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P29762)

Gene Name	CRABP1
About Gene Symbol	CRABP1 (Cellular Retinoic Acid Binding Protein 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

