

KD-Validated Anti-Janus Kinase 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1086

Product Overview

Name	KD-Validated Anti-Janus kinase 1 Rabbit Monoclonal Antibody
Accession(Primary)	P23458
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human JAK1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 133 kDa , observed, 133 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P23458)

Gene Name	JAK1
About Gene Symbol	JAK1 (Janus Kinase 1) is a receptor tyrosine kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

