

KD-Validated Anti-Membrane Metalloendopeptidase Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1099

Product Overview

Name	KD-Validated Anti-Membrane Metalloendopeptidase Rabbit Monoclonal Antibody
Accession(Primary)	P08473
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CD10
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 86 kDa, observed, 100 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000


Target information (P08473)

Gene Name	MME
About Gene Symbol	MME (Membrane Metalloendopeptidase) is a protease enzyme involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-CD10 antibody (Cat#ABG1099). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-CD10 antibody (Cat#ABG1099, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-CD10 antibody (Cat#ABG1099). CD10 expression in wild type (WT) and CD10 shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-CD10 antibody (Cat#ABG1099, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.