

KD-Validated Anti-Transcription Elongation factor A1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1101

Product Overview


Name	KD-Validated Anti-Transcription elongation factor A1 Rabbit Monoclonal Antibody
Accession(Primary)	P23193
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human TCEA1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 34 kDa, observed, 34 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P23193)


Gene Name	TCEA1
About Gene Symbol	TCEA1 (Transcription Elongation Factor A1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HepG2 cells with Transcription elongation factor A1 antibody (Cat#ABG1101, 1:1,000). Nuclei were stained blue with DAPI; Transcription elongation factor A1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 μ m.