

KD-Validated Anti-STUB1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1103


Product Overview

Name	KD-Validated Anti-STUB1 Rabbit Monoclonal Antibody
Accession(Primary)	Q9UNE7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human STUB1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 35 kDa , observed, 35 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (Q9UNE7)

Gene Name	STUB1
About Gene Symbol	STUB1 (STIP1 Homology And U-Box Containing Protein 1) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway and adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and





Immunocytochemical staining of C2C12 cells with STUB1 antibody (Cat#ABG1103, 1:1,000). Nuclei were stained blue with DAPI; STUB1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.