

KD-Validated Anti-Carbamoyl-phosphate synthase 1 Rabbit

Rabbit monoclonal antibody


ABG1107

Product Overview



Name	KD-Validated Anti-Carbamoyl-phosphate synthase 1 Rabbit Monoclonal Antibody
Accession(Primary)	P31327
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CPS1
Species Reactivity	Human
Calculated MW	Predicted, 165 kDa , observed, 150 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (P31327)


Gene Name	CPS1
About Gene Symbol	CPS1 (Carbamoyl-Phosphate Synthase 1) is a biosynthetic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and ident



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC