

KD-Validated Anti-Adenylosuccinate Lyase Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG1110

Product Overview


Name	KD-Validated Anti-Adenylosuccinate Lyase Mouse Monoclonal Antibody
Accession(Primary)	P30566
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human ADSL
Species Reactivity	Human, Rat
Calculated MW	Predicted, 55 kDa, observed, 48 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,500-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P30566)


Gene Name	ADSL
About Gene Symbol	ADSL (Adenylosuccinate Lyase) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC