

KD-Validated Anti-TXNDC5 Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG1114

Product Overview


Name	KD-Validated Anti-TXNDC5 Mouse Monoclonal Antibody
Accession(Primary)	Q8NBS9
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human TXNDC5
Species Reactivity	Human
Calculated MW	Predicted, 48 kDa, observed, 48 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,500-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q8NBS9)


Gene Name	TXNDC5
About Gene Symbol	TXNDC5 (Thioredoxin Domain Containing 5) is a protein kinase and signaling regulator involved in angiogenesis. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechan



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC

Thioredoxin domain containing 5-Alexa Fluor® 647