

KD-Validated Anti-F-Box Protein 11 Mouse Monoclonal

Mouse monoclonal antibody


ABG1116

Product Overview

Name	KD-Validated Anti-F-Box Protein 11 Mouse Monoclonal Antibody
Accession(Primary)	Q86XK2
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human FBXO11
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 104 kDa, observed, 120 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,500-1:5,000 Flow Cytometry (FCM): 1:2,000


Target information (Q86XK2)

Gene Name	FBXO11
About Gene Symbol	FBXO11 (F-Box Protein 11) is a methyltransferase involved in ubiquitin-proteasome pathway. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-FBXO11 antibody (Cat#ABG1116). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-FBXO11 antibody (Cat#ABG1116, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-FBXO11 antibody (Cat#ABG1116). FBXO11 expression in wild-type (WT) and FBXO11 knockdown (KD) HeLa cells with 30 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-FBXO11 antibody (Cat#ABG1116, 1:1,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of FBXO11 expression in HeLa cells using anti-FBXO11 antibody (Cat#ABG1116, 1:2,000). Green, isotype control; red, FBXO11.