

KD Validated Anti-Phospho-HSF1 (S326) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1119

Product Overview


Name	KD Validated Anti-Phospho-HSF1 (S326) Rabbit Monoclonal Antibody
Accession(Primary)	Q00613
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-HSF1 (S326)
Species Reactivity	Human
Calculated MW	Predicted, 57 kDa, observed, 82 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q00613)


Gene Name	HSF1
About Gene Symbol	HSF1 (Heat Shock Transcription Factor 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC