

KD-Validated Anti-PPP3CA Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1124

Product Overview

Name	KD-Validated Anti-PPP3CA Rabbit Monoclonal Antibody
Accession(Primary)	Q08209
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Calcineurin A
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 59 kDa , observed, 59 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (Q08209)

Gene Name	PPP3CA
About Gene Symbol	PPP3CA (Protein Phosphatase 3 Catalytic Subunit Alpha) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mec

