

KD-Validated Anti-Complement C3 Rabbit Monoclonal

Rabbit monoclonal antibody

ABG1135

Product Overview

Name	KD-Validated Anti-Complement C3 Rabbit Monoclonal Antibody
Accession(Primary)	P01024
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human C3 / C3a
Species Reactivity	Human
Calculated MW	Predicted, 187 kDa , observed, 187 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (P01024)

Gene Name	C3
About Gene Symbol	C3 (Complement C3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

