

KD-Validated Anti-CC2D1A Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1147

Product Overview


Name	KD-Validated Anti-CC2D1A Rabbit Monoclonal Antibody
Accession(Primary)	Q6P1N0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from CC2D1A
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 104 kDa , observed, 130 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q6P1N0)


Gene Name	CC2D1A
About Gene Symbol	CC2D1A (Coiled-Coil And C2 Domain Containing 1A) is a protein kinase and signaling regulator involved in NF-κB signaling and vesicular trafficking. It is an important research target in cancer, cell signaling, and therapeutic target research and is widel



Western blotting analysis using anti-CC2D1A antibody (Cat#ABG1147). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-CC2D1A antibody (Cat#ABG1147, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-CC2D1A antibody (Cat#ABG1147). CC2D1A expression in wild type (WT) and CC2D1A shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-CC2D1A antibody (Cat#ABG1147, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Flow cytometric analysis of CC2D1A expression in C2C12 cells using CC2D1A antibody (Cat#ABG1147, 1:2,000). Green, isotype control; red, CC2D1A.