

KD-Validated Anti-CIRBP Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1148


Product Overview

Name	KD-Validated Anti-CIRBP Rabbit Monoclonal Antibody
Accession(Primary)	Q14011
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from CIRBP
Species Reactivity	Human
Calculated MW	Predicted, 19 kDa , observed, 19 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q14011)

Gene Name	CIRBP
About Gene Symbol	CIRBP (Cold Inducible RNA Binding Protein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with CIRBP antibody (Cat#ABG1148, 1:1,000). Nuclei were stained blue with DAPI; CIRBP was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 μ m.