

KD-Validated Anti-ACSL4 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1150

Product Overview


Name	KD-Validated Anti-ACSL4 Rabbit Monoclonal Antibody
Accession(Primary)	O60488
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from ACSL4
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 79 kDa, observed, 73 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (O60488)


Gene Name	ACSL4
About Gene Symbol	ACSL4 (Acyl-CoA Synthetase Long Chain Family Member 4) is a biosynthetic enzyme involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blotting analysis using anti-ACSL4 antibody (Cat#ABG1150). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-ACSL4 antibody (Cat#ABG1150, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-ACSL4 antibody (Cat#ABG1150). ACSL4 expression in wild type (WT) and ACSL4 shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-ACSL4 antibody (Cat#ABG1150, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Flow cytometric analysis of ACSL4 expression in HT-1080 cells using ACSL4 antibody (Cat#ABG1150, 1:2,000). Green, isotype control; red, ACSL4.



Immunocytochemical staining of HT-1080 cells with ACSL4 antibody (Cat#ABG1150, 1:1,000). Nuclei were stained blue with DAPI; ACSL4 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.