

KD-Validated Anti-CLIP1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1157


Product Overview

Name	KD-Validated Anti-CLIP1 Rabbit Monoclonal Antibody
Accession(Primary)	P30622
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from CLIP170
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 162 kDa , observed, 170 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P30622)

Gene Name	CLIP1
About Gene Symbol	CLIP1 (CAP-Gly Domain Containing Linker Protein 1) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with CLIP1 antibody (Cat#ABG1157, 1:1,000). Nuclei were stained blue with DAPI; CLIP1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.