

KD-Validated Anti-Caveolin 3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1159


Product Overview

Name	KD-Validated Anti-Caveolin 3 Rabbit Monoclonal Antibody
Accession(Primary)	P56539
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from Caveolin-3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 17 kDa, observed, 17 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P56539)

Gene Name	CAV3
About Gene Symbol	CAV3 (Caveolin 3) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HT-1080 cells with Caveolin 3 antibody (Cat#ABG1159, 1:1,000). Nuclei were stained blue with DAPI; Caveolin 3 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.