

KD-Validated Anti-Cathepsin D Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1171

Product Overview


Name	KD-Validated Anti-Cathepsin D Rabbit Monoclonal Antibody
Accession(Primary)	P07339
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Cathepsin D
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 45 kDa , observed, 28, 46 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P07339)

Gene Name	CTSD
About Gene Symbol	CTSD (Cathepsin D) is a lysosomal cathepsin protease involved in vesicular trafficking. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC