

KD-Validated Anti-IKK α Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1176

Product Overview

Name	KD-Validated Anti-IKK α Rabbit Monoclonal Antibody
Accession(Primary)	O15111
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human IKK α
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 85 kDa , observed, 85 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (O15111)

Gene Name	CHUK
About Gene Symbol	CHUK (Component Of Inhibitor Of Nuclear Factor Kappa B Kinase Complex) is a protein kinase and signaling regulator involved in NF- κ B signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely st

