

KD-Validated Anti-Caspase 3 p12 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1180


Product Overview

Name	KD-Validated Anti-Caspase 3 p12 Rabbit Monoclonal Antibody
Accession(Primary)	P42574
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Caspase-3 p12
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 32 kDa , observed, 12,29,33 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P42574)

Gene Name	CASP3
About Gene Symbol	CASP3 (Caspase 3) is a cysteine protease (caspase) involved in apoptosis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of Caspase 3 p12 expression in C2C12 cells using Caspase 3 p12 antibody (Cat#ABG1180, 1:2,000). Green, isotype control; red, Caspase 3 p12.