

KD-Validated Anti-CD99 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1188

Product Overview

Name	KD-Validated Anti-CD99 Rabbit Monoclonal Antibody
Accession(Primary)	P14209
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CD99
Species Reactivity	Human
Calculated MW	Predicted, 19 kDa , observed, 28 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P14209)

Gene Name	CD99
About Gene Symbol	CD99 (CD99 Molecule (Xg Blood Group)) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

