

KD-Validated Anti-ADP Ribosylation factor 6 Rabbit Monoclonal Antibody


ABG1190

Product Overview


Name	KD-Validated Anti-ADP ribosylation factor 6 Rabbit Monoclonal Antibody
Accession(Primary)	P62330
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human ARF6
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 20 kDa, observed, 18 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P62330)


Gene Name	ARF6
About Gene Symbol	ARF6 (ADP Ribosylation Factor 6) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HT-1080 cells with ADP ribosylation factor 6 antibody (Cat#ABG1190, 1:1,000). Nuclei were stained blue with DAPI; ADP ribosylation factor 6 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.