

# KD-Validated Anti-Chloride intracellular channel 4 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1203

## Product Overview

|                             |                                                                                                              |
|-----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                 | KD-Validated Anti-Chloride intracellular channel 4<br>Rabbit Monoclonal Antibody                             |
| <b>Accession(Primary)</b>   | Q9Y696                                                                                                       |
| <b>Host Species</b>         | Rabbit                                                                                                       |
| <b>Clonality</b>            | Monoclonal                                                                                                   |
| <b>Isotype</b>              | Rabbit IgG                                                                                                   |
| <b>Immunogen</b>            | A synthesized peptide derived from human CLIC4                                                               |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                                                            |
| <b>Calculated MW</b>        | Predicted, 29 kDa , observed, 29 kDa                                                                         |
| <b>Storage Buffer</b>       | Supplied in PBS (pH 7.4) containing 50% glycerol,<br>and 0.02% sodium azide.                                 |
| <b>Recommended Dilution</b> | Western Blotting (WB): 1:1,000-1:5,000 Flow<br>Cytometry (FCM): 1:2,000 Immunocytochemistry<br>(IC): 1:1,000 |


## Target information (Q9Y696)

|                          |                                                                                                                                                                                                                                                                                       |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CLIC4                                                                                                                                                                                                                                                                                 |
| <b>About Gene Symbol</b> | CLIC4 (Chloride Intracellular Channel 4) is a membrane transporter involved in membrane transport. It is an important research target in pharmacology, drug resistance, and membrane biology and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |




Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Chloride intracellular channel 4 antibody (Cat#ABG1203). Total cell lysates (30  $\mu$ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Chloride intracellular channel 4 antibody (Cat#ABG1203, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.




Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Chloride intracellular channel 4 antibody (Cat#ABG1203). Chloride intracellular channel 4 expression in wild type (WT) and Chloride intracellular channel 4 shRNA knockdown (KD) HeLa cells with 30  $\mu$ g of total cell lysates.  $\beta$ -Tubulin serves as a loading control. The blot was incubated with anti-Chloride intracellular channel 4 antibody (Cat#ABG1203, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of Chloride intracellular channel 4 expression in HAP-1 cells using Chloride intracellular channel 4 antibody (Cat#ABG1203, 1:2,000). Green, isotype control; red, Chloride intracellular channel 4.



Immunocytochemical staining of HAP-1 cells with Chloride intracellular channel 4 antibody (Cat#ABG1203, 1:1,000). Nuclei were stained blue with DAPI; Chloride intracellular channel 4 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20  $\mu$ m.