

KD-Validated Anti-CSDE1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1206


Product Overview

Name	KD-Validated Anti-CSDE1 Rabbit Monoclonal Antibody
Accession(Primary)	O75534
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CSDE1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 89 kDa , observed, 89 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (O75534)

Gene Name	CSDE1
About Gene Symbol	CSDE1 (Cold Shock Domain Containing E1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of C2C12 cells with CSDE1 antibody (Cat#ABG1206, 1:1,000). Nuclei were stained blue with DAPI; CSDE1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.