

KD-Validated Anti-ABHD5 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1209

Product Overview

Name	KD-Validated Anti-ABHD5 Rabbit Monoclonal Antibody
Accession(Primary)	Q8WTS1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human ABHD5
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 39 kDa , observed, 45 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q8WTS1)

Gene Name	ABHD5
About Gene Symbol	ABHD5 is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

