

KD-Validated Anti-BCAR1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1218

Product Overview

Name	KD-Validated Anti-BCAR1 Rabbit Monoclonal Antibody
Accession(Primary)	P56945
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human BCAR1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 93 kDa , observed, 130 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P56945)

Gene Name	BCAR1
About Gene Symbol	BCAR1 (BCAR1 Scaffold Protein, Cas Family Member) is a signaling scaffold protein involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

