

KD-Validated Anti-Casein Kinase 2 beta Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1219

Product Overview


Name	KD-Validated Anti-Casein kinase 2 beta Rabbit Monoclonal Antibody
Accession(Primary)	P67870
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Casein Kinase 2 beta
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 25 kDa , observed, 25 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 immunocytochemistry (IC): 1:100-1:1,000

Target information (P67870)


Gene Name	CSNK2B
About Gene Symbol	CSNK2B (Casein Kinase 2 Beta) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Geniu Biotechnologies LLC



Copyright ©2024 Geniu Biotechnologies LLC



Copyright ©2024 Geniu Biotechnologies LLC