

KD-Validated Anti-CHAF1B Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1222

Product Overview

Name	KD-Validated Anti-CHAF1B Rabbit Monoclonal Antibody
Accession(Primary)	Q13112
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human p60 CAF1
Species Reactivity	Human
Calculated MW	Predicted, 61 kDa , observed, 70 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (Q13112)

Gene Name	CHAF1B
About Gene Symbol	CHAF1B (Chromatin Assembly Factor 1 Subunit B) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

