

KD-Validated Anti-Core-binding factor subunit beta Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1223

Product Overview

Name	KD-Validated Anti-Core-binding factor subunit beta Rabbit Monoclonal Antibody
Accession(Primary)	Q13951
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human CBFb
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 22 kDa , observed, 22 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000


Target information (Q13951)

Gene Name	CBFB
About Gene Symbol	CBFB (Core-Binding Factor Subunit Beta) is a transcription factor involved in gene transcription and regulation. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Core-binding factor subunit beta antibody (Cat#ABG1223). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Core-binding factor subunit beta antibody (Cat#ABG1223, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Core-binding factor subunit beta antibody (Cat#ABG1223). Core-binding factor subunit beta expression in wild type (WT) and Core-binding factor subunit beta shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-Core-binding factor subunit beta antibody (Cat#ABG1223, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of Core-binding factor subunit beta expression in HAP-1 cells using Core-binding factor subunit beta antibody (Cat#ABG1223, 1:2,000). Green, isotype control; red, Core-binding factor subunit beta.