

KD-Validated Anti-Menin 1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1236

Product Overview

Name	KD-Validated Anti-Menin 1 Rabbit Monoclonal Antibody
Accession(Primary)	O00255
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Menin
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 68 kDa; observed, 68 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western blotting (WB): 1:1,000-1:5,000Flow Cytometry (FCM): 1:2,000


Target information (O00255)

Gene Name	MEN1
About Gene Symbol	MEN1 (Menin 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-menin 1 antibody (Cat#ABG1236). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-menin 1 antibody (Cat#ABG1236, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-menin 1 antibody (Cat#ABG1236). Menin 1 expression in wild type (WT) and menin 1 (MEN1) knockdown (KD) HSHC cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-menin 1 antibody (Cat#ABG1236, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of Menin 1 expression in HepG2 cells using anti-Menin 1 antibody (Cat#ABG1236, 1:2,000). Green, isotype control; red, Menin 1.