

KD-Validated Anti-Helicase like transcription factor Rabbit

Rabbit monoclonal antibody


ABG1239

Product Overview

Name	KD-Validated Anti-Helicase like transcription factor Rabbit Monoclonal Antibody
Accession(Primary)	Q14527
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human HLTF
Species Reactivity	Human
Calculated MW	Predicted, 114 kDa , observed, 120 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000


Target information (Q14527)

Gene Name	HLTF
About Gene Symbol	HLTF (Helicase Like Transcription Factor) is a transcription factor involved in chromatin remodeling and epigenetics and ubiquitin-proteasome pathway. It is an important research target in gene regulation, development, and cancer research and is widely st




Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Helicase like transcription factor antibody (Cat#ABG1239). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Helicase like transcription factor antibody (Cat#ABG1239, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-Helicase like transcription factor antibody (Cat#ABG1239). Helicase like transcription factor expression in wild type (WT) and helicase like transcription factor shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-Helicase like transcription factor antibody (Cat#ABG1239, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Immunocytochemical staining of HepG2 cells with Helicase like transcription factor antibody (Cat#ABG1239, 1:1,000). Nuclei were stained blue with DAPI; Helicase like transcription factor was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: High. Scale bar: 20 µm.