

KD-Validated Anti-LMAN1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1247

Product Overview


Name	KD-Validated Anti-LMAN1 Rabbit Monoclonal Antibody
Accession(Primary)	P49257
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human LMAN1
Species Reactivity	Human, Mouse, Rat
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (P49257)


Gene Name	LMAN1
About Gene Symbol	LMAN1 (Lectin, Mannose Binding 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blotting analysis using anti-LMAN1 antibody (Cat#ABG1247). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-LMAN1 antibody (Cat#ABG1247, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-LMAN1 antibody (Cat#ABG1247). LMAN1 expression in wild type (WT) and LMAN1 shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-LMAN1 antibody (Cat#ABG1247, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Flow cytometric analysis of LMAN1 expression in HepG2 cells using LMAN1 antibody (Cat#ABG1247, 1:2,000). Green, isotype control; red, LMAN1.



Immunocytochemical staining of HepG2 cells with LMAN1 antibody (Cat#ABG1247, 1:1,000). Nuclei were stained blue with DAPI; LMAN1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.