

KD-Validated Anti-MCM3 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1256

Product Overview

Name	KD-Validated Anti-MCM3 Rabbit Monoclonal Antibody
Accession(Primary)	P25205
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human MCM3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 91 kDa , observed, 100 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P25205)

Gene Name	MCM3
About Gene Symbol	MCM3 (Minichromosome Maintenance Complex Component 3) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic

