

KD-Validated Anti-Matrix metalloproteinase 14 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1260

Product Overview

Name	KD-Validated Anti-Matrix metalloproteinase 14 Rabbit Monoclonal Antibody
Accession(Primary)	P50281
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human MMP14
Species Reactivity	Human, Rat
Calculated MW	Predicted, 66 kDa , observed, 65 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P50281)

Gene Name	MMP14
About Gene Symbol	MMP14 (Matrix Metalloproteinase 14) is a protease enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

