

KD-Validated Anti-Phospho-RSK1(S380) Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG1281

Product Overview


Name	KD-Validated Anti-Phospho-RSK1(S380) Rabbit Monoclonal Antibody
Accession(Primary)	Q15418
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Phospho-RSK1(S380)
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 83 kDa, observed, 74 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:1,000

Target information (Q15418)


Gene Name	RPS6KA1
About Gene Symbol	RPS6KA1 (Ribosomal Protein S6 Kinase A1) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied




Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Immunocytochemical staining of HeLa cells with Phospho-RSK1(S380) antibody (Cat#ABG1281, 1:1,000). Nuclei were stained blue with DAPI; Phospho-RSK1(S380) was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.